



Simax Media Bottles

SIMAX®

Made from heat resistant borosilicate 3.3 glass in accordance with ISO 3585 and conforming to ISO 4796-1. Complete with blue plastic cap and pouring ring resistant to 140°C. GL45 thread. Autoclavable.

Code	Alt Ref	Capacity, mL	Cap Size	Pack
BOT53020	2070/50	50	GL32	10
BOT53040	2070/100	100	GL45	10
BOT53060	2070/250	250	GL45	10
BOT53080	2070/500	500	GL45	10
BOT53100	2070/1000	1000	GL45	10
BOT53120	2070/2000	2000	GL45	10
BOT5314	2070/5000	5000	GL45	Each
BOT5316	2070/10000	10,000	GL45	Each
BOT5318	2070/20000	20,000	GL45	Each



DURAN Polypropylene Caps and Rings

DURAN®

Caps for all DURAN® glass screw thread bottles with GL threads. Features moulded internal cone/lip seal - no need to separate for cleaning and sterilisation. Made from durable polypropylene.

Polypropylene pouring cap rings provide a sharp edge to the DURAN bottle neck allowing a clean pour and avoiding messy and hazardous dribbles. Not required for DURAN bottles with GL25 thread sizes.

Suitable for automatic cleaning, chemical sterilisation and autoclaving. Temperature range -40°C to +140°C. Autoclavable up to 140°C, steam heat. Unsuitable for DURAN soda-lime bottles. Please visit www.scientificlabs.co.uk for alternative colour caps.

Supplied in packs of 10

Code	Alt Ref	Description	Cap size	Colour
BOT5230	2923913	Cap	GL25	Blue
BOT5232	2923919	Cap	GL32	Blue
BOT5234	2923928	Cap	GL45	Blue
293382802	293382802	Cap	GL45	Yellow
29338286	29338286	Cap	GL45	Green
29338288	29338288	Cap	GL45	Grey
BOT5236	2924219	Pouring ring	GL32	Blue
BOT5238	292422809	Pouring ring	GL45	Blue
1089911	1089911	Pouring ring	GL45	Green
1089914	1089914	Pouring ring	GL45	Grey
1089917	1089917	Pouring ring	GL45	Yellow