

Hach Photometric Tests Cont.

Code	Alt Ref	Type	Measuring Range	Method
WAT3649	LCK338	Total Nitrogen, LatoN®	20-100 TN b	2,6-Dimethyl phenol
WAT3640	LCK049	ortho - Phosphorus (inc ortho-Phosphate)	1.6-30mg/L PO ₄ -5-90mg/L PO ₄	Vanadate-molybdate
WAT3653	LCK348	Phosphate total (inc ortho-Phosphate)	0.5-5 PO ₄ -P, 1.5-15 PO ₄	Phosphomolybdenum blue
WAT3655	LCK350	Phosphate total (inc ortho-Phosphate)	2-20 PO ₄ -P, 6-60 PO ₄	Phosphomolybdenum blue
WAT3657	LCK383	TOC	5-50 TOC	Expulsion method
WAT3665	LCK384	TOC	50-500 TOC	Expulsion method
WAT3667	LCK380	TOC	2-65 TOC	Difference method

Lovibond® Comparator 2000+

The Lovibond® Comparator system is an extremely versatile, modular system. Compact and portable, it is simple to use yet uncompromising in terms of precision. The glass colour standards are fade-free, resistant to chemicals, scratchproof and compliant to international standards, eg ISO 7393/2.

- Colorimetric, visual water analysis
- Versatile, modular system
- 400 different test discs
- Precision results
- Optional lighting unit
- Accuracy of Lovibond® Reagents



2
YEARS

Code	Alt Ref	Description
WAT6000	142000	Comparator 2000+
WAT6016	142050	Daylight Unit light source
WAT2098	171010	Comparator daylight 2000+ light source
WAT6300	230070	Ammonia 3/113 disc
WAT6302	233930	Chlorine 3/40/K disc
WAT6304	230270	Phosphate 3/133 disc
WAT6306	230690	Manganese 3/169 disc
WAT2080	234010	Chlorine disc 3/40A
WAT2082	221130	Phenol Red disc 2/1J
WAT2084	234060	Chlorine disc 3/40E
WAT2086	237000	Phosphate disc 3/71
WAT2088	232070	Chlorine disc 3/2ARP
WAT2096	511751BT	Phenol red tablets
WAT2058	411120	AF112A free and combined Cl kit
WAT2060	411090	AF110 high and low chlorine Kit
WAT2062	413570	AF357 drinking water kit
WAT2164	145505	Checkit 10mL plastic cells
WAT2166	354243	Glass moulded cell 10mL + cap