

## Tristel Medistel High Level Instrument Disinfectant

Tristel

Medistel is a disinfectant for flexible and rigid endoscopes and surgical instrument based on quaternary ammonium compounds. It is primarily designed for open sink use but is also suitable for use in ultrasonic baths. Medistel is bactericidal, fungicidal, virucidal and tuberculocidal. Demonstrated efficacy against Gram- and Gram+ bacteria, yeasts, fungi, viruses, mycobacteria and spores according to European Norms.

Code	Alt Ref	Description	Pack	Price
TRI1375	MD311	Concentrate bottle with dosing cap	6 x 1L	£209.00



## Rely+On™ Virkon™ Disinfectant

The broad spectrum formulation of Rely+On™ Virkon™ is unique. No other disinfectant has the same powerful composition or extensive portfolio of performance and safety testing data. It combines application flexibility with broad spectrum efficacy, on hard surfaces, and in the face of organic challenge. These qualities make Rely+On™ Virkon™ the disinfectant of choice for use in medical facilities, pathology and biosafety containment laboratories.

- Independently proven highly effective against over 100 strains of virus in 22 viral families, over 400 strains of bacteria, and over 60 strains of fungi and yeast using a wide variety of contact times, temperatures and organic challenge levels
- Unrivalled spectrum
- Built-in colour indicator of strength
- Resistant to inactivation by organic material
- Contains a surfactant - combines cleaning and disinfection in one operation
- Avoids the problems of pre-cleaning which spreads infection
- Tablets are convenient to store and easy to handle



Code	Pack	Makes	Price
CLE1550	50 x 50g sachets	1 sachet = 5L disinfectant	£177.00
CLE1552	6 x 500g	1 pack = 50L disinfectant	£180.00
CLE1554	1 x 5kg	1 drum = 500L disinfectant	£180.00
33012	20 x 10 x 5gm tablets	10 x 5gm tablets = 5L disinfectant	£105.00
CLE1556	12 x 50 x 5gm tablets	50 x 5gm tablets = 25L disinfectant	£337.00

## Klercide 70/30 Sterile Alcohol

70% denatured, pharma grade ethanol blended with 30% deionised water. Supplied sterile by gamma irradiation and double bagged. Cleanroom processed and 0.2µm filtered. Bactericidal and fungicidal to EN 1650, EN 1278 and EN13697.

Code	Pack	Price
3035340	6 x 1L	£158.00



## Sterillium® Disinfectant

Alcohol based hand disinfectant with skin protection factor in easy dispense bottle. It is effective against a broad range of microorganisms and viruses, and fast acting with residual effect.

- Reduces transient germs by 99.99% within 30 seconds
- Also effective against skin flora
- Especially mild to the skin because of skin care additives
- Clinical, microbial, virucidal, toxicological and cosmetic

Code	Pack	Price
CLE1580	500mL	£16.72
CLE1582	1L	£18.61

