

## Pyrex® Cylinder Heavy Duty Squat Form Class B

PYREX®

- Complies with ISO 4788
- Manufactured from Pyrex® borosilicate glass with robust hexagonal foot
- Graduations and inscriptions in blue enamel
- Heavy duty construction for increased mechanical strength
- Squat form for convenience, especially in confined spaces
- Heavily beaded rim and spout

Supplied in packs of 2

Code	Alt Ref	Capacity, mL	Graduations	Tolerance
CYL13200	3285-10	10	1	0.3
CYL13220	3285-25	25	1	0.5
CYL13240	3285-50	50	1	1
CYL13260	3285-100	100	2	1
CYL13280	3285-250	250	5	2
CYL13300	3285-500	500	10	5
CYL13320	3285-1L	1000	20	10



## SLS SELECT Cylinder Class B



With spout, glass foot and white enamel graduations. Hexagonal base. To DIN12680 (ISO4788).

Code	Capacity, mL
CYL1070	5
CYL1072	10
CYL1074	25
CYL1076	50
CYL1078	100
CYL1080	250
CYL1082	500
CYL1084	1000



## Glass Cylinder Stoppered

Fitted with interchangeable plastic stopper. Manufactured in accordance with BS EN ISO 4788:2005.

Code	Capacity, mL	Graduations, mL	Stopper
CYL1100	5	0.1	10/13
CYL1102	10	0.2	10/13
CYL1104	25	0.5	10/13
CYL1106	50	1.0	14/15
CYL1108	100	1.0	14/15
CYL1110	250	2.0	19/17
CYL1112	500	5.0	24/20
CYL1114	1000	10.0	34/35
CYL1116	2000	20.0	34/35

