

## Buchi E-812 and E-816 Hot Extraction Systems

Programmable, electronically controlled unit for solid-liquid extraction for fat determination following the hot extraction principle (e.g. Randall). Suitable for extractions with petroleum ether, diethyl ether, hexane or chloroform. Available as 2-place [E-812 HE] or 6-place unit [E-816 HE].


- Automated hot extraction according to Randall or Goldfisch
- High solvent recover due to in built cooled tank
- Safe - Cooling water, level of solvent and heat source are directly monitored, the complete apparatus is specially produced to view complete tightness of joints while a protective shield provides an extra precaution
- Easy operation - fully programmable



Code	Alt Ref	Description
EVA4290	049100	Buchi E-812 Extraction unit for hot extraction, two place
EVA4292	047580	Buchi E-816 Extraction unit for hot extraction, six place

## Buchi E-812 HE and E-816 HE selected accessories

Code	Alt Ref	Description
EVA4234	034925	Quartz sand for extraction and hydrolysis 1 x 25 kg
EVA4236	037689	Quartz sand for extraction and hydrolysis 1 x 2.5 kg
EVA4296	045674	Set of 12 beakers for Buchi HE
EVA4298	047643	6 position Beaker holder for Buchi equipment
EVA4300	049426	Set of two beakers for Buchi HE
EVA4250	049430	Set of 2 x glass sample tubes with frit for Buchi equipment
EVA4252	051903	Holder for 6 x glass sample tubes, PTFE
EVA4254	11055738	Set of 12 x glass sample tubes with frit for Buchi equipment



## Can't find what you're looking for?

Let us know! We have over 300,000 products available, so if you can't find it our website will. Visit today.

[www.scientificlabs.co.uk](http://www.scientificlabs.co.uk)

