



Gilson PIPETMAN® G Pipette

PIPETMAN G incorporates a completely redesigned and lubricated piston system.

- Extremely low pipetting force
- Extensive range covers volume requirements from 0.2µL up to 10mL
- PIPETMAN Comfort Handle can clip on for enhanced ergonomics and personalisation

PIPETMAN® G Starter Kit

- Three PIPETMAN G pipettes: P20G, P200G, and P1000G (fitted with stainless steel ejectors)
- Two racks of certified-quality PIPETMAN Diamond® Tips: D200 and D1000
- Three adhesive-backed SINGLE™ pipette hangers
- One 'Gilson Guide to Pipetting' Handbook
- 3 interchangeable plastic tip ejectors

PIPETMAN® G Micro Volume Kit

- Three PIPETMAN G pipettes: P2G, P10G, and P100G (fitted with stainless steel ejectors) Dual tip ejector adapters for P2G and P10G
- Two racks of Diamond tips DL10 and D200
- Three adhesive-backed SINGLE™ pipette hangers
- One 'Gilson Guide to Pipetting' Handbook
- 3 interchangeable plastic tip ejectors

PIPETMAN® G 4 Pipette Kit

- 4 PIPETMAN G Pipettes: P2G, P20G, P200G, and P1000G
- 3 racks of certified-quality PIPETMAN DIAMOND Tips: D10L, D200 and D1000
- SINGLE™ pipette holders
- 1 Gilson Guide to Pipetting

Code	Description	Volume Range, µL	Ejector
F144054P	P2G	0.2-2.0	Plastic
F144055P	P10G	0.5-10	Plastic
F144056P	P20G	2-20	Plastic
F144057P	P100G	10-100	Plastic
F144058P	P200G	20-200	Plastic
F144059P	P1000G	100-1000	Plastic
F144054M	P2G	0.2-2.0	Metal
F144055M	P10G	0.5-10	Metal
F144056M	P20G	2-20	Metal
F144057M	P100G	10-100	Metal
F144058M	P200G	20-200	Metal
F144059M	P1000G	100-1000	Metal
F144066	P5000G	500-5000	-
F144067	P10mLG	1000-10,000	-
F167900	PIPETMAN® G Starter Kit	-	-
F167800	PIPETMAN® G Microvolume Kit	-	-
F167360	PIPETMAN® G 4 pipette kit	-	-



CARROUSEL, Stands and Reservoirs



Suitable for use with the PIPETMAN P, PIPETMAN L, PIPETMAN G, PIPETMAN Neo and Microman.

Code	Description
F161401	CARROUSEL™ Pipette Stands - Holds up to 7 Pipettes
F161406	SINGLE™ Pipette Holder
F161405	TRIO™ Pipette Stand - Holds up to 3 pipettes
F267660	Reagent reservoir 25mL - 10 packs of 10 reservoirs
F267670	Reagent reservoir 50mL - 10 packs of 10 reservoirs