


Volac Class AS Bulb Pipettes, One Mark

These class AS Bulb pipettes satisfy International Standards (BS ISO 648) for accuracy.

Supplied in packs of 5

Code	Capacity, mL	Tolerance, \pm mL	Works Certificate	Price
PIP12000	1	0.008	No	£59.00
PIP12300	1	0.008	Yes	£42.00
PIP12020	2	0.01	No	£41.00
PIP12320	2	0.01	Yes	£42.00
PIP12040	3	0.015	No	£24.00
PIP12340	3	0.015	Yes	£44.00
PIP12060	4	0.015	No	£36.00
PIP12360	4	0.015	Yes	£50.00
PIP12080	5	0.015	No	£47.00
PIP12380	5	0.015	Yes	£46.00
PIP12100	10	0.02	No	£38.00
PIP12400	10	0.02	Yes	£52.00
PIP12120	15	0.025	No	£42.00
PIP12140	20	0.03	No	£44.00
PIP12440	20	0.03	Yes	£58.00
PIP12160	25	0.03	No	£45.00
PIP12460	25	0.05	Yes	£59.00
PIP12180	50	0.05	No	£60.00
PIP12480	50	0.08	Yes	£74.00
PIP1220	100	0.08	No	£42.00
PIP1250	100	0.08	Yes	£58.00



MBL® Class B Pipettes

- Complies with ISO 835
- Graduations in permanent amber stain
- Manufactured from soda-lime glass with robust toolled jets

MBL®



Type 1

Calibrated to deliver from zero at the top to any graduation line down to the shoulder.

Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Graduations, mL	Tolerance, \pm mL	Colour Code	Price
PIP22000	PRA760B	1	0.01	0.01	Yellow	£12.44
PIP22020	PRA762B	2	0.02	0.02	Black	£14.11
PIP22040	PRA764B	5	0.05	0.05	Red	£14.55
PIP22060	PRA766B	10	0.10	0.10	Orange	£15.33
PIP22080	PRA768B	25	0.10	0.20	Green	£22.50