

SLS Select Enviroscreen Environmental and Hygiene Sampling Swabs Continued



Surface Sampling: Premoistened Swab

- Ready to use and cost effective sampling swabs.
- Perfect for hygiene monitoring of surfaces
- Consistent sampling with no mess
- Surfaces should be wiped post sampling
- Sterilised by gamma irradiation
- 2 year shelf life

Code	Description	Pack	Price
SWA1002	Premoistened in neutralising buffer, viscose tip, polystyrene shaft	250	£78.00
SWA1004	Premoistened in neutralising buffer, viscose tip, blue B/P PS shaft	250	£79.00
SWA1006	Premoistened in deionised water, viscose tip, polystyrene shaft, tube, peel pouch	100	£27.80



Meat Carcass Sampling: Sponge Kits

- Ready-to-use all in one kit for the environmental sampling of meat carcasses
- No additional preparation required
- Can be used for the recovery of total viable count (TVC) and Enterobacteriaceae (including *Salmonellae* and *E. coli*)
- Designed according to Food Standard Agency recommendations and complies with EU regulation no. 2073/2005.
- Easy to use, reliable and cost effective
- 50cm² blue cellulose sponge
- Pre-moistened with choice of solutions
- Flexible breakpoint handle available - Envirostik™
- Comply with EU regulation no. 2073/2005
- Sterilised by gamma irradiation
- 1 year shelf life



Code	Description	Pack	Price
SWA1038	50cm ² blue sponge, 0.9% NaCl, 100mL plastic container	95	£252.00
SWA1143	50cm ² blue sponge, dosed in saline (0.9% NaCl) easy open stomacher pouch	200	£240.00
SWA1145	50cm ² blue sponge, dosed in peptone saline diluent (MRD), easy open stomacher pouch	200	£223.00
SWA1149	50cm ² blue sponge, dosed in buffered peptone water, easy open stomacher pouch	200	£234.00

Surface Sampling: MRD System (No Neutralisers)

- Ideal for non-disinfected areas.
- Sensitive enough for hand sampling
- Quick turnaround analysis, 4 to 6 hours
- Sterilised by gamma irradiation
- 2 year shelf life

Code	Description	Pack	Price
SWA1124	Liquid MRD medium, viscose tip, polystyrene shaft, tube in peel pouch	100	£81.00

