

## Hamilton 1700 Series Microlitre Syringes

HAMILTON

The 1700 series syringes are the gastight version of the original 700 series syringes. Other 1700 series syringes and accessories are available, please enquire for further information.

### Fixed Needle, Sharp



Code	Type	Capacity, $\mu\text{L}$	Price
SYR7090	1701	10	£44.50
SYR7092	1702	25	£53.00
SYR7050	1705	50	£57.00
SYR7200	1725	250	£60.00
SYR7202	1750	500	£59.00


### Removable Needle, Sharp



Code	Type	Capacity, $\mu\text{L}$	Price
SYR7094	1701	10	£63.00
SYR7095	1702	25	£68.00
SYR7060	1705	50	£68.00
SYR7204	1725	250	£99.00
SYR7206	1750	500	£73.00

### Luer Tip

Code	Type	Capacity, $\mu\text{L}$	Price
SYR7096	1701	10	£73.00
SYR7098	1702	25	£62.00
SYR7070	1705	50	£64.00
SYR7072	1710	100	£62.00
SYR7208	1725	250	£66.00
SYR7074	1750	500	£63.00



## Hamilton 7000 Series Microlitre Syringes

HAMILTON

The 7000 series syringes are designed on the dead-volume-free principal of positive displacement. A fine tungsten plunger extends through the whole length of the syringe to the needle end.

### Fixed Needle, Sharp

Code	Type	Capacity, $\mu\text{L}$	Price
SYR7220	7101	1	£83.00
SYR7222	7102	2	£110.00
SYR7224	7105	5	£98.00



### Fixed Needle, Blunt

Code	Type	Capacity, $\mu\text{L}$	Price
SYR7226	7101	1	£86.00
SYR7228	7102	2	£112.00
SYR7230	7105	5	£101.00

