

# Eppendorf ThermoMixer C



## FREE SmartBlock

Eppendorf ThermoMixer® C combines unique mixing performance with excellent temperature control to guarantee complete, dependable and reproducible test results. Improves user assay results by mixing and incubating samples at the same time. With quick and easy-to-exchange SmartBlocks for all common lab vessels, new and improved program functions, and a high level of stability.

- ✓ Mixing Frequency: 300 - 3000 RPM
- ✓ Temperature Range: 15°C below ambient to 100°C
- ✓ Temperature Accuracy: ±0.5 °C at 20 – 45°C
- ✓ Dims, w x d x h, mm: 206 x 304 x 136

Select one of the following FREE SmartBlocks when you purchase the ThermoMixer C (Shipped Separately)

Code	For tubes
F9392	8 x 5.0mL
F9394	24 x 1.5mL
F9396	24 x 0.5mL
F9398	24 x 2.0mL
F9400	4 x 50mL
F9402	8 x 15mL
F9404	24 Cryotubes 1.5-2mL

Code	Description	List Price, £	List Price, €
F9390	ThermoMixer C without Block	2097.00	2621.25

Worth €543.75

