



MICROMAN E

Ergonomic Pipetting For Problem Liquids

Guarantee accuracy, precision and comfort when pipetting volatile and viscous samples. MICROMAN E together with disposable capillary pistons ensures contamination-free pipetting and protects your pipette from corrosion.

50%
OFF

- Experience complete and rapid pipetting of viscous, dense, hot or cold samples.
- Avoid leaks when pipetting volatile (high vapour pressure) liquids.
- No cross-contamination or residual sample carry-over.

MICROMAN E		
MODEL VOLUME	ORDER REF	OFFER PRICE
M10E 1-10 µL	FD10001PR18	£248.00 £124.00
M25E 3-25 µL	FD10002PR18	
M50E 20-50 µL	FD10003PR18	
M100E 10-100 µL	FD10004PR18	
M250E 50-250 µL	FD10005PR18	
M1000E 100-1000 µL	FD10006PR18	

Not to be used in conjunction with any other offer or contract discount.
Valid until 31/12/18 to UK customers only.

How To Order:

Call 0115 982 1111 or 01236 431857 for Scotland

or visit www.scientificlabs.co.uk

or email slsinfo@scientific-labs.com





CAPILLARY PISTONS (CPs)

100% Protection Against Contamination

MICROMAN Capillary Pistons (CPs) offer 100% protection against contamination. Available as standard CPs or sterile CPs which are free of DNA, RNA, DNAses, RNAses, ATP, pyrogens and trace metals. Now with 35% off list price.



**35%
OFF CPs**

- Large volume range
1µL to 1000 µL across six different models.
- Free of dyes, no risk of contamination from CP components.
- Lightweight and slim CPs that fit in even the narrowest and deepest tubes.

How To Order:

Call 0115 982 1111 or 01236 431857 for Scotland

or visit www.scientificlabs.co.uk

or email slsinfo@scientific-labs.com

CAPILLARY PISTONS - 35% OFF				
CAPILLARY PISTON	ORDER REF	QTY/BOX	LIST PRICE	OFFER PRICE
CP10	F148412PR18	192	£60.00	£39.00
	F148312PR18	960	£237.00	£154.05
CP25	F148112PR18	200	£50.00	£32.50
	F148012PR18	960	£261.00	£169.65
CP50	F148113PR18	200	£50.00	£32.50
	F148013PR18	960	£261.00	£169.65
CP100	F148414PR18	192	£60.00	£39.00
	F148314PR18	960	£237.00	£154.05
CP250	F148114PR18	200	£50.00	£32.50
	F148014PR18	960	£261.00	£169.65
CP1000	F148560PR18	182	£60.00	£39.00
STERILISED CAPILLARY PISTONS - 35% OFF				
CP10ST	F148413PR18	192	£68.00	£44.20
	F148313PR18	960	£337.00	£219.05
CP25ST	F148712PR18	576	£202.00	£131.30
CP50ST	F148713PR18	576	£202.00	£131.30
CP100ST	F148415PR18	192	£68.00	£44.20
	F148315PR18	960	£337.00	£219.05
CP250ST	F148714PR18	576	£202.00	£131.30
CP1000ST	F148180PR18	192	£68.00	£44.20

Not to be used in conjunction with any other offer or contract discount.
Valid until 31/12/18 to UK customers only.